

Background

Cryo-EM aids characterization of sensitive battery interfaces

Cryogenic electron microscopy (cryo-EM) is a **valuable tool for characterizing sensitive materials**, as the low temperatures minimize many electron beam-induced artifacts & damage. In lithium-ion battery research, this has enabled transmission & scanning transmission electron microscopy (TEM & STEM) characterization of electrode-electrolyte interphases at high spatial resolution.

Cryo-EM sample prep is a challenge; samples must be:

- (1) **Frozen quickly** to preserve native structures
- (2) **Thin** enough for S/TEM characterization
- (3) **Protected from air/moisture** (for batteries).

Traditional sample prep methods are ex situ and leave significant time between sample cycling & freezing:




Fig. 1: Typical cryo-EM sample prep workflow for batteries.

This **limits time resolution**, allowing structural relaxation, diffusion, or other processes to occur between cycling & cryo-EM, and **fails to capture interfacial structures in their electrified state**.

Goal: Time-resolved method to capture in-situ structures at electrified interfaces.

Methods

In-situ freezing of electrified state




Fig. 2: Freezing samples in-situ under applied bias (or current) allows cryo-EM of structures correlating with times A, B, & C.

1. Integrating plunge freezing with in-situ electrochemical testing

Sample holder developed for Thermo Fisher Vitrobot connects to current/voltage source, allows sample to be plunged into liquid nitrogen (LN₂) at programmed points in electrochemical testing.

$\Delta t_s \approx 1.5$ s, for step-based plunging (e.g. "plunge at 1.5 V")

$\Delta t_t \approx 0.1$ s, for time-based plunging (e.g. "plunge after 10 s")



Fig. 3: Biasing sample holder in Vitrobot.

2. Modified coin cell for cold retrieval of active material




Fig. 4: (a) Retrieving sample, (b) modified cap & cell.

Electrode, on TEM grid, retrieved while continuously submerged in LN₂, for **fast, protected transfer**.

Results

Test case: Silicon nanoparticle half cell with GenF electrolyte

Half cells: Si nanoparticle (NP)-coated TEM grid or Si wafer grid working electrode, Li metal counter electrode, Celgard separator, & GenF electrolyte ([1.2 M LiPF₆ in 3:7 w/w EC:EMC] + 3 wt.% FEC).

Control samples show:

- (1) **Consistent electrochemistry** between modified or conventional cells, cycled in in-situ holder or conventional cyler
- (2) **Suitability for multimodal cryo-EM**, with vitrified electrolyte around active materials.




Fig. 5: (a) Cryo-electron diffraction shows electrolyte is vitrified. (b) Si NPs within frozen electrolyte. (c) Cryo-energy dispersive X-ray spectroscopy of native oxide on Si wafer electrode.

Impact: Preservation of electrified state for cryo-EM characterization

By integrating plunge freezing with electrochemical testing of battery materials, we have achieved **3-4 orders of magnitude improvement in time resolution** over conventional cryo-EM sample prep techniques and **preservation of interfaces in their electrified state**.

Future work can **extend this method to other systems** with sensitive materials & dynamic structure under bias, such as perovskite photovoltaics or biohybrid catalysts.

Acknowledgements: The authors acknowledge Sang-Don Han for insight into modified coin cell designs and Gerard Michael Carroll for preparation of Si nanoparticles. The authors also acknowledge Steve Robbins, Rory Andrykowski, and Jeff Alleman, as well as John Watt and the Center for Integrated Nanotechnologies, for support with Vitrobot modifications.